

Introduction to TOPAS

Profile fitting and Rietveld analysis

Holger Cordes, Ph.D.

Applications Scientist, XRD

5465 E. Cheryl Parkway

Madison, WI 53711-5373

608-276-3048

Holger.Cordes@bruker-axs.com

Outline:

TOPAS Introduction

TOPAS Interface

Profile fitting

- Fitting with analytical functions
- Convolution based approaches and Fundamental Parameters
- Empirical determination of instrument functions
- Whole Powder Pattern Decomposition
- Obtaining information about micro structure
- Anisotropic Refinement models
- Determining Degree of crystallinity with Peaks phases
- Determining Degree of crystallinity with hkl phases

- Quantitative Analysis with Rietveld
 - Introduction
 - Corrections
 - Examples:
 - Round robin example with preferred orientation
 - Cement sample
 - Mixture of minerals, ...
 - PONKCS

- Generalized software for profile and structure analysis

Seamless integration of all currently employed profile fit techniques and related applications

- Single Line Fitting
- Indexing
(LSI, LP-Search)
- Whole Powder Pattern Decomposition
(Pawley, Le Bail)
- Ab initio structure determination
(Simulated annealing, charge flipping, Fourier analysis)
- Structure refinement
(Rietveld)
- Quantitative phase analysis
(Rietveld, PONKCS)

Size / Strain
Degree of crystallinity

TOPAS Interface:


GUI and Launch mode

■ TOPAS has a graphical user-interface (GUI) and a text based (Launch)mode


- GUI supports most Profile fitting and Rietveld applications
- Beginners should use this mode of operation.

■ LAUNCH MODE

- More features for advanced users
- For structure solution
- For automation or batch operation



TOPAS GUI: User Interface



Parameter window
with instrument and
structure information

Fit window

wt.% for each phase

Measured curve
Calculated curve
Difference curve

TOPAS GUI: Parameter Window



Import raw data files

Data files will be displayed
in treeview.

Instrument settings have
to be entered, if Fundamental
Parameters (peak type FP) is
used

Add a peak phase (for profile
fitting), an hkl-phase (for
lattice parameter refinements)
or load a structure (for
quantitative Rietveld
refinements)

Parameter values

Toggle any parameter between Fix and
Refine in „Code“ column

	Use	Value	Code	Error	Min	Max
Background						
Chebyshev	<input checked="" type="checkbox"/>		@			
Order		1				
1/X Bkg	<input type="checkbox"/>	1000	Refine	0		
Goniometer radii						
Primary radius (mm)		217.5				
Secondary radius (mm)		217.5				
Equatorial Convolutions						
Point detector	<input checked="" type="checkbox"/>					
Receiving Slit Width (mm)	<input checked="" type="checkbox"/>	0.1	Fix	0		
FDS Shape, angle(°)	<input checked="" type="checkbox"/>	0.5	Fix	0		
Beam spill, sample length	<input type="checkbox"/>	50	Fix	0		
VDS irradiated length (mm)	<input type="checkbox"/>	12	Fix	0		
VDS Scale Intensity	<input type="checkbox"/>					
Capillary	<input type="checkbox"/>					
Linear PSD	<input type="checkbox"/>					
Tube Tails	<input type="checkbox"/>					
Axial Convolutions						
Full Axial Model	<input checked="" type="checkbox"/>					
Source length (mm)		12	Fix	0		
Sample length (mm)		15	Fix	0		
RS length (mm)		12	Fix	0		
Prim. Soller (°)	<input checked="" type="checkbox"/>	2.3	Fix	0		
Sec. Soller (°)	<input checked="" type="checkbox"/>	2.3	Fix	0		
N Beta		30				
Finger_et_al	<input type="checkbox"/>					
Simple Axial Model (mm)	<input type="checkbox"/>	12	Fix	0		

TOPAS GUI:

Fit window




Start Refinement (Step by step)

Refinement until convergence

Residual of the Refinement at every refinement step

graphical output of the residual as a functions of iteration number

Launch mode on/off



Profile Fitting Methods

How do you describe a Peak?

A peak may be described by the following peak parameters:

- background (Due to air scatter, fluorescence, overlapping peak tails, diffuse scatter)
- Position (integral intensity depends more on weight fraction, atomic position or other relevant physical parameters than intensity maximum)
- Intensity
- Line width (FWHM)
- Shape (How fast does the intensity drop off?)
- Asymmetry (FWHM, Left and Right)

Profile Fitting Methods


Analytical approach

■ Analytical Profile Fitting

- Characterized by peak shape functions (PSFs) with an explicit and relatively simple mathematical form which can be differentiated **analytically** with respect to each of the refined parameters
- Most common are the Gaussian, Lorentzian, sums of Gaussians or Lorentzians Voigt, pseudo Voigt (PV) and PearsonVII functions (PVII)
- Convenient to use
- Large number of parameters required to fit profiles well over the whole 2θ range pattern frequently leading to
 - Correlation problems
 - Loss of uniqueness
 - Instability of the refinement procedure

Analytical approach

Example: CeO₂



- Import the raw file and zoom on the main low angle peak
- Open Peak Details window (F3)
- Choose SPVII (Split Pearson VII) peak type and drop a peak on the main peak with the mouse. This adds a Peaks Phase
- Start refinement (F6)
- Get peak parameters

File: C:\Topas 4.2\tutorial\CeO2\CeO2.raw

Profile Fitting Methods

Convolution based methods

■ Direct Convolution Approach

- Discrimination of instrument and sample contributions:

$$Y(2\theta) = (W \times G) \times S, \text{ where } I = (W \times G)$$

W: Emission profile **G**: Instrument contribution **S**: Sample contribution

TOPAS offers 3 direct convolution approaches:

1. **Empirical parametrization** of line profile shapes by convolution of appropriate functions
 - Profile parameters have NO physical meaning
2. Explicit discrimination of instrument and sample contributions:
Measured Instrument Function
 - Sample related profile parameters have physical meaning: Size/strain/...
3. Explicit discrimination of instrument and sample contributions:
Calculated Instrument Function
 - All profile parameters have a physical meaning
 - **Fundamental Parameters Approach (Peak Type FP)**

Convolution approach

Empirical parametrization




- Convolute any appropriate functions to achieve a best fit
 - Minimize the number of functions / function parameters to
 - minimize parameter correlation
- An excellent approach for all profile fit applications using any instrument, if micro-structure information is NOT of interest (Refined profile parameters have NO physical meaning)
- Knowledge of the most common contributions to line profile shapes and their dependence on angle helps:

<i>Contribution</i>	<i>Convolution</i>	<i>Angular dependence</i>
Detector (slit)	Hat	Constant
Crystallite size	Lorentzian	$1/\cos(\Theta)$
Strain	Gaussian	$\tan(\Theta)$
Axial divergence	Circles	$-1/\tan(\Theta)$

Convolution approach

Empirical parametrization Examples



Size-Strain Round Robin, Balzar (2004); "sharp data".


Convolution approach

Exercise: Empirical parametrization



File: C:\Topas 4.2\tutorial\misc\Lab6_a1.raw
Emission profile: CuK α 2.lam

- Use Peak type FP for any convolution based approach
- For the empirical approach make sure that no instrument parameters are used (uncheck all)
- Add Peaks Phase, uncheck crysize and select "Additional Convolutions"
- "Add a Convolution"
- Select a predefined Convolution function, define a 2Th dependency and set Value to refine
- Add several convolutions until you achieve a good fit
- Trial and Error approach



Convolution approach

Measured instrument function approach

Step 1: Determine an instrument function

- Obtain data from a suitable standard:
 - Reflection: LaB6 (SRM 660a)
 - Transmission: Si (SRM 640a)
- Convolute any appropriate functions to achieve a best fit (Empirical parameterization!)
- The number of functions / function parameters is irrelevant, the more the better!
- Fix all refineable profile parameters and save an instrument file

Step 2: Refine the actual data using the instrument function

- Load the instrument file
- Refine on micro-structure parameters as required

Convolution approach

Measured instrument function


File: C:\Topas 4.2\examples\GADDS LaB6 long scan.raw

Emission profile: CuK α 2.lam

Determine instrument function

Determination of Gadds instrument parameters using "LaB6 long scan.raw"

hkl_SG 221 a=4.15 0.00 %




Convolution approach

Measured instrument function

	Use	Value	Code	Error
Goniometer radii				
Primary radius (mm)		250		
Secondary radius (mm)		250		
Equatorial Convolutions				
Point detector	<input checked="" type="checkbox"/>			
Receiving Slit Width (mm)	<input type="checkbox"/>	7.103858e-(Refine		0
FDS Shape, angle(°)	<input checked="" type="checkbox"/>	1.195284	Refine	0
Beam spill, sample length (mm)	<input type="checkbox"/>	50	Fix	0
VDS irradiated length (mm)	<input type="checkbox"/>	10	Fix	0
VDS Scale Intensity	<input type="checkbox"/>			
Capillary	<input type="checkbox"/>			
Linear PSD	<input type="checkbox"/>			

	Conv. Type	2Th Dependence	Use	Value	Code	Error	Min
1	Lorentzian	1/Cos(Th)	<input checked="" type="checkbox"/>	0.01972032	Refine	0	
2	Circles	-1/Tan(Th)	<input checked="" type="checkbox"/>	0.01379071	Refine	0	

- Any "slit width" or "function" under the Divergent beam tab or any additional convolution can be refined to get the best possible fit for the instrument peak shape.
- No limit to number of functions or parameters, the more the better the fit.
- After refining LaB₆ these parameters they are set to "fixed" and used as an empirical instrument function for refinement of unknown samples.




Use an "unknown" sample with identical instrument settings
File: C:\Topas 4.2\examples\GADDS Anatase 5pc Rutile.raw


Convolution approach

Calculated Instrument Function


Fundamental Parameters




W



G



S



$$Y(2\theta) = (W \times G) \times S$$

Convolution approach

Calculated Instrument Function

Fundamental Parameters



Implemented only for Bragg-Brentano Geometry

Notation for slit settings:


Instrument parameter notation in TOPAS in Divergent beam tab:

RS: Receiving Slit (entered in mm)

FDS: Fixed Divergence Slit (entered in degrees)


VDS Length: Variable Slit (mm correspond to beam width on sample)

Full Axial Model: enter primary and secondary Soller slit sizes here



Single Line Fitting

Choose an Emission Profile



TOPAS - [Lab6_a1.raw]

File View Fit Launch Tools Window Help

Global
Lab6_a1.raw
Emission Profile
Background
Instrument
Corrections
Miscellaneous

Load Emission Profile
Save Emission Profile
Add Emission Line
Paste INP to Node/Selections


	Area	WL (Å)	Lortz. HW (nGauss)	HW (Å)
1	0.653817	1.540596	0.501844	0
2	0.346183	1.544493	0.626579	0

	Use	Value
Ymin on Ymax		0.001
No Th dependence	<input type="checkbox"/>	
For LAM cursor	<input type="checkbox"/>	
Lam for Bragg angle	<input type="checkbox"/>	0
Calculate Lam	<input type="checkbox"/>	

- Predefined emission profiles for most common sealed tube systems and radiations are predefined in the directory `c:\topas3\lam\`
- The default profile `Cuka2.lam` has 2 lines, $\text{Cu}\alpha_1$ and $\text{Cu}\alpha_2$.
- If a residual of the $k\beta$ line is visible for measurements with Ni filter and very intense peaks, an emission line with 1.39222\AA can be added and its intensity be refined.
- For data with very sharp peaks the profile `Cuka5.lam` will give slightly better refinements.
- Check the „For LAM cursor“ box to visualize the emission profile with the cursor
- For optics that may distort the wavelength spectrum (e.g. Goebel mirrors, Polycaps), the emission profile can be refined.

File: C:\Topas 4.2\tutorial\misc\Lab6_1.raw

Modeling the background in TOPAS



TOPAS - [Lab6_a1.raw]

File View Fit Launch Tools Window Help

Global

Lab6_a1.raw

Emission Profile

Background

Instrument

Corrections

Miscellaneous

Background Rpt/Text

	Use	Value	Code	Error
Chebychev	<input checked="" type="checkbox"/>		@	
Order		3		
1/X Bkg	<input checked="" type="checkbox"/>	1000	Refine	0


Chebychev Coefficients

	Coefficient	Error
1	0	0
2	0	0
3	0	0
4	0	0

Paste INP to Node/Selections

- Most backgrounds can be fitted with a polynomial function (Chebychev Polynomial) and relatively few parameters
- The 1/x function is used to fit the rise the background at low angles caused by scatter from the direct beam
- Too few (<3) or too many parameters (>6) can be a source for serious errors, especially if broad peaks are present
- If an amorphous halo is present, it can often be fitted by entering a peak phase with a single peak

Entering instrument parameters in TOPAS



- Highlight the instrument item in the treeview below the scan file
- In the tab „divergent beam“ the slit sizes for Bragg-Brentano geometry are entered in either mm or degrees.
- Radius depends on measurement circle of diffractometer
 - D4, D500, D5000: 200.5mm
 - D8: usually 217.5 or 250mm
- For Point detector: Check the use column for point detector and enter settings for receiving and divergence slit
- For Linear PSDs (LynxEye and Vantec-1) check here and enter angular 2θ range and divergence
- Soller slits sizes are entered under Full Axial model
 - Default size for D8: 2.5°
 - Default size for D4: 4°
- Instrument settings can be saved, loaded or made a default setting

FPA Example for LaB₆ (111)




Single line fitting: Exercise



W

$W \times G_{Eq} \times G_{Ax}$

$(W \times G_{Eq} \times G_{Ax}) \times S$



File: C:\Topas 4.1\tutorial\misc\Lab6_a1.raw

Emission profile: CuKα2.lam

Background order 3

Radius: 217 mm

Divergence: 0.5°

Receiving slit: 0.1mm


Soller slits: 2.3°

Advantages of FPA:

- Identical instrument contribution to every peak
- Fewer parameters: Faster, more stable computation
- No parameter turn-on sequence for refinements
- Average crystallite size/strain for each phase

Profile Fitting in TOPAS

Peaks Phases



Adding a Peaks Phase will add a phase without any position or intensity constraints

This is typically used to determine accurate peak positions for Indexing, peak areas or micro-structural parameters

Select the Peak Details window (F3)

Several Peak Types are available

FP: Fundamental Parameters

PV: Pseudo Voigt function

PVII: Pearson VII function

SPV: Split Pseudo Voigt function

SPVII: Split Pearson VII function


The second column allows display of position, area, crystallite size or strain at the peak

Select the peak type SPVII and click with the mouse near the peak maxima of your scan

File: C:\Topas 4.1\tutorial\misc\quartz.raw

Profile Fitting in TOPAS

Peaks Phases, Peak types



Select SPVII peak type and click near the maximum of the three peaks with the mouse.


Run the Refinement

Select individual peak display

The third peak has almost 50% more intensity than the others, even though the maximum intensity is visually lower

Also note the considerably higher asymmetry

This is caused by parameter correlation and can be avoided by using constraints that ensure the same peak width for all peaks



Values		Codes	Errors	Min				
Max	Scherrer	LVoHB (nm)	LVoF-PWHM (nm)	e0	Rpt/Text			
Use	Position	Area	FWHM Left	FWHM Right	M Left	M Right		
1	✓ 67.73133	511.6396	0.05504582	0.03754237	2.506625	1.278807		
2	✓ 68.13123	586.8861	0.04567504	0.03597622	2.915417	9.179479		
3	✓ 68.29791	707.9003	0.06775916	0.04041264	0.7614877	1.232083		

File: C:\Topas 4.1\tutorial\misc\quartz.raw

Emission profile: CuKα₂_analyt.lam

Background order 3

Profile Fitting in TOPAS

Peaks Phases, Constraining parameters



Values Codes Errors Min

Max Scherrer LVolHB (nm) LVol-FWHM (nm) e0 Rpt/Text

	Use	Position	Area	FWHM Left	FWHM Right	M Left	M Right
1	<input checked="" type="checkbox"/>	67.73133	511.6396	0.05	0.05	1	1
2	<input checked="" type="checkbox"/>	68.13123	586.8861	0.05	0.05	1	1
3	<input checked="" type="checkbox"/>	68.29791	707.9003	0.05	0.05	1	1

Values Codes Errors Min

Max Scherrer LVolHB (nm) LVol-FWHM (nm) e0 Rpt/Text

	Use	Position	Position	Area	FWHM Left	FWHM Right	M Left	M Right
1	<input checked="" type="checkbox"/>	67.73133	Refine	Refine	FWHML	FWHMR	ML	Mr
2	<input checked="" type="checkbox"/>	68.13123	Refine	Refine	FWHML	FWHMR	ML	Mr
3	<input checked="" type="checkbox"/>	68.29791	Refine	Refine	FWHML	FWHMR	ML	Mr


To change the value of a whole column to the same, click on the upper right corner of the grid

Type identical start value for FWHM and M

Select Codes tab and enter the parameter names for FWHM and M

Run the Refinement again

The peaks now have identical shapes and correct intensities



Values Codes Errors Min

Max Scherrer LVolHB (nm) LVol-FWHM (nm) e0 Rpt/Text

	Use	Position	Area	FWHM Left	FWHM Right	M Left	M Right
1	<input checked="" type="checkbox"/>	67.72707	532.3152	0.04967582	0.04142547	1.86854	1.421101
2	<input checked="" type="checkbox"/>	68.13088	768.883	0.04967582	0.04142547	1.86854	1.421101
3	<input checked="" type="checkbox"/>	68.29765	385.4827	0.04967582	0.04142547	1.86854	1.421101

File: C:\Topas 4.1\tutorial\misc\quartz.raw

Emission profile: CuK α 2_analyt.lam

Background order 3

Profile Fitting in TOPAS

Peaks Phases, Fundamental parameters as an intrinsic constraint



TOPAS - C:\users\8043 files for topas training\test.pro - [Quartz.raw]

File View Fit Launch Tools Window Help

Global
Quartz.raw
Emission Profile
Background
Instrument
Corrections
Miscellaneous
Peaks Phase:0
Peaks FP

Divergent beam Additional Convolutions Rpt/Text

	Use	Value
Goniometer radii		
Primary radius (mm)		217.5
Secondary radius (mm)		217.5
Equatorial Convolutions		
Point detector	<input checked="" type="checkbox"/>	
Receiving Slit Width (mm)	<input checked="" type="checkbox"/>	0.1
FDS Shape, angle(°)	<input checked="" type="checkbox"/>	1
Beam spill, sample length (mm)	<input type="checkbox"/>	50
VDS irradiated length (mm)	<input type="checkbox"/>	12
VDS Scale Intensity	<input type="checkbox"/>	
Capillary	<input type="checkbox"/>	
Linear PSD	<input type="checkbox"/>	
Tube Tails	<input type="checkbox"/>	
Axial Convolutions		
Full Axial Model	<input checked="" type="checkbox"/>	

Global
Quartz.raw
Emission Profile
Background
Instrument
Corrections
Miscellaneous
Peaks Phase:0
Peaks FP

Values Codes Errors Min Max LVolHB (nm) LVolFW (nm)

	Type	Use	Position
1	FP	<input checked="" type="checkbox"/>	67.7318
2	FP	<input checked="" type="checkbox"/>	68.13618
3	FP	<input checked="" type="checkbox"/>	68.30136

Change the Peak type to FP
Enter Instrument parameters for that scan.

Instrumental contribution to all peaks is identical
→ intrinsic constraint, however anisotropic effects can be refined individually

Difference between instrumental width and measured width is caused by microstructure such as size broadening or strain

File: C:\Topas 4.1\tutorial\misc\quartz.raw

Emission profile: CuKa2_analyt.lam

Background order 3

Radius: 217.5mm

Receiving slit 0.1mm

FDS 1°

4° Soller slits

Global
Quartz.raw
Emission Profile
Background
Instrument
Corrections
Miscellaneous
Peaks Phase:0
Peaks FP

Values Codes Errors Min Max LVolHB (nm) LVolFWHM (nm) e0 Rpt/Text

	Use	Position	Area	Cry size L (nm)	Cry size G (rStrain L)	Strain G
1	<input checked="" type="checkbox"/>	67.7318	595.1021	237.0	200.0	0.1
2	<input checked="" type="checkbox"/>	68.13618	833.056	285.8	200.0	0.1
3	<input checked="" type="checkbox"/>	68.30136	444.1307	228.4	200.0	0.1

Profile Fitting in TOPAS

Peaks Phases

TOPAS - C:\users\8043 files for topas training\test.pro - [Quartz.raw]

File View Fit Launch Tools Window Help

Quartz.raw

- Emission Profile
- Background
- Instrument
- Corrections
- Miscellaneous
- Peaks Phase:0
- Peaks SPVII

Values Codes Errors Min Max

LVolHB (nm) LVol-FWHM (nm) e0 Additional Convolutions

	Type	Use	Position	Area	Cry size L (n	Cry size G
1	SPVII	<input checked="" type="checkbox"/>	67.72927	532.1825	200.0	200.0
2	SPVII	<input checked="" type="checkbox"/>	68.13253	769.1925	200.0	200.0
3	SPVII	<input checked="" type="checkbox"/>	68.29935	386.5809	200.0	200.0

Hide peak sticks

Save Phase

Save Peaks as DIF file

Change to d_{hkl} phase

Create Indexing Range

Create Pawley range

Delete Peaks Phase

Paste INP to Node/Selections

Hide peak sticks in display

Save Peaks as an input file or dif file

Create Indexing Range for Indexing single phase samples

Creates an duplicate range with identical instrument parameters but hkl phase



File: C:\Topas 4.2\tutorial\misc\quartz.raw

Emission profile: CuKa2_analyt.lam

Background order 3

Profile Fitting in TOPAS

Peaks Phases, Calculating Areas for crystallinity determination



File: C:\Topas

4.2\tutorial\doc\m2a.raw

Emission profile: CuKα2_analyt.lam


Background order 1

1/x function

- Add two peak phases, rename them amorphous and crystalline
- Use one SPV peak for the amorphous hump
- Select crystalline and insert a peak for each reflection (Verify that you are not fitting the background with these peaks otherwise use symmetric functions or constraints)
- Highlight Miscellaneous and Calculate Degree of Crystallinity

Profile Fitting in TOPAS

hkl Phases for Lattice Parameter Refinement



File: C:\Topas 4.2\tutorial\ce02\ce02.raw

Measurement circle radius: 173mm

Emission profile: CuKα5.lam

Background order 3, 1/x function



Divergence: 1°, Receiving slit: 0.1mm

Soller slits: prim. 5.1°, sec. 8.6°

- Add an hkl phase
- Enter a space group number or symbol or pick a group from the space group editor
- Enter approx. lattice parameter (5.4 for CeO₂) and set to refine
- Enter instrument parameters
- Check zero error correction under Corrections
- Run refinement

Profile Fitting in TOPAS

hkl Phases for Lattice Parameter Refinement

	h	k	l	m	d	th2	Use	I	Code	Error
1	1	1	1	8	3.12353	28.55412	✓	706	Refine	1.79
2	0	0	2	6	2.70506	33.08923	✓	259	Refine	1.26
3	0	2	2	12	1.91277	47.49592	✓	993	Refine	3.37
4	3	1	1	24	1.63121	56.35754	✓	1.05e+003	Refine	3.99
5	2	2	2	8	1.56177	59.10516	✓	215	Refine	1.91
6	0	0	4	6	1.35253	69.43404	✓	254	Refine	2.34
7	3	3	1	24	1.24117	76.72346	✓	658	Refine	3.95
8	0	4	2	24	1.20974	79.09934	✓	447	Refine	3.33
9	4	2	2	24	1.10434	88.45695	✓	713	Refine	4.43
10	5	1	1	24	1.04118	95.43439	✓	337	Refine	2.51e+008

- Refined lattice parameter
- Refined average crystallite size (Scherrer equation)
- Peak list with intensities under the hkl Is tab

Profile Fitting Methods in TOPAS

Exercise: Refine Lattice parameter of Ibuprofen



File: C:\Topas
4.1\tutorial\Examples\ibuprofen\ibuprofen.raw
 Emission profile: CuK α 1.lam (Vario Data)
 1mm Capillary data with LynxEye detector

Obtain space group and start values for lattice parameters from EVA

EVA - [JCP2.2CA:00-032-1723]

File Pattern View Options Window Help

PDF Number: 00-032-1723 View Quality: High


C₁₃H₁₈O₂

Ibuprofen

General Comments Authors Additional Subfiles


Cell Parameters		Crystal Data	
Latt.: Monoclinic		Molecular weight: 206.28	
S.G.: P2 ₁ /c (14)		Volume [CD]: 1226.33	
a: 14.667	alpha: 99.46	Dx: 1.117	
b: 7.893	beta: 99.46	Dm:	
c: 10.731	gamma: 4	I/Cor:	
a/b: 1.85682	Z: 4		
c/b: 1.35853			

Phase Details					
Additional Convolutions					
Rpt/Text					
	Use	Value	Code	Error	
Use Phase	<input checked="" type="checkbox"/>				
Le Bail	<input type="checkbox"/>				
Delete hkl on Refinement	<input checked="" type="checkbox"/>				
LP Search	<input type="checkbox"/>	0.4			
Spacegroup		p2 ₁ /c			
a (Å)		14.6571964	Refine	0.0002981	
b (Å)		7.8834464	Refine	0.0001460	
c (Å)		10.7215003	Refine	0.0001447	
beta (°)		99.42143	Refine	0.00082759	



Peak Shape corrections

Absorption in low density or loosely packed samples



- Identical Ibuprofen samples measured with identical scan parameters but different sample holders
- Flat low background holder
- Low background holders with 0.4mm deep cavity
- PMMA holder with 1mm depth
- For low density samples absorption correction is necessary

Corrections			
Cylindrical sample (Sabine)			
Rpt/Text			
	Use	Value	Code
Peak shift			
Zero error	<input type="checkbox"/>	0	Refine
Sample displacement (mm)	<input type="checkbox"/>	0	Refine
Intensity Corrections			
LP factor	<input type="checkbox"/>	0	Fix
Surface Rghnss Pitschke e	<input type="checkbox"/>		
Surface Rghnss Suortti	<input type="checkbox"/>		
Sample Convolutions			
Absorption (1/cm)	<input checked="" type="checkbox"/>	100	Refine
Sample Thickness (mm)	<input type="checkbox"/>	10	Refine
Scale Intensity	<input type="checkbox"/>		
Sample Tilt (mm)	<input type="checkbox"/>	0	Refine

File: C:\Topas 4.2\examples\Ibuprofen D2.raw

Emission profile: CuKa2.lam

Try out Absorption correction

Profile Fitting Methods in TOPAS

Exercise: Determine degree of crystallinity in Ibuprofen testmixture



File: C:\Topas 4.2\Examples\testmixture D2.raw and ...\Ibuprofen D2.raw

Emission profile: CuK α 2.lam


D2 Phaser with 141mm radius

2.5° primary Soller and 4° secondary soller

Lynxeye with 2.5° opening

Both scans were collected with identical conditions

Obtain degree of crystallinity of testmixture D2.raw (assume that all excipients are amorphous and Ibuprofen D2.raw scan is 100% crystalline)



Microstructural information

Size and Strain



TOPAS default (in GUI)

Size/Strain Component functions (Double Voight Approach)

■CrySize_L: $\text{lor_fwhm} = 0.1 \lambda / \text{cs1} \cos(\text{Th});$

■CrySize_G: $\text{gauss_fwhm} = 0.1 \lambda / \text{cs2} \cos(\text{Th});$

■Strain_L: $\text{lor_fwhm} = \text{str1} \tan(\text{Th});$

■Strain_G: $\text{gauss_fwhm} = \text{str2} \tan(\text{Th});$

cs1, cs2, str1, str2: Refineable parameters, arbitrary parameter names!




Microstructural information

Size/Strain



File: C:\Topas 4.2\tutorial\misc\Lab6_1.raw

Lvol: Volume weighted mean column height
Size is shape dependent: Simple relationship for cubes and spheres but not anything else



Microstructural information

Size/Strain



Common assumption:

- Lorentzian part is Size broadening and Gaussian part is Strain.

Better:

- Model each broadening mechanism with Voight functions (Balzar and Ledbetter 1995)
- If the data is measured in large range, the component can still be separated by angular dependencies ($\tan\Theta$ for strain and $1/\cos\Theta$ for size)
- If any Gaussian size broadening is present, it will lead to better fit
- Size distribution can be attempted in Launch Mode by calculating the integral breadth of Lorentzian (Cauchy) size contribution. Works in simple cases and if instrumental function is well known.

General Considerations

- In most applications, the objective is primarily to improve the fit, with an expectation to improve indexing, structure determination and structure refinement
- Anisotropic micro-structure parameters can be derived
- Example functions
 - Spherical harmonics (Popa, 1998)
 - 2nd rank tensors (Le Bail & Jouanneaux, 1997)
 - Multi-dimensional distribution of lattice metrics (Stephens, 1999)

■ TOPAS Tutorial Examples:

- Norbonane Fitch & Jovic, 1993
- LT-ZrMo₂O₈ Allen et al., 2003
- NaBenzo R. Dinnebier et al., 1999
- Clay S. Hillier, 2003
- SerpMix1 M. Raudsepp, 2004

File: C:\Topas 4.2\tutorial\anisoLS\...

■ Keywords that can be hkl dependent

- *th2_offset* E (anisotropic peak shifts - macrostrain)
- *scale_pks* E (preferred orientation)

- *lor_fwhm* E (anisotropic line width)
- *gauss_fwhm* E (anisotropic line width)
- *hat* E (anisotropic line width)
- *one_on_x_conv* E (anisotropic line asymmetry)
- *exp_conv_const* E (anisotropic line asymmetry)
- *circles_conv* E (anisotropic line asymmetry)
- *user_defined_convolution*

■ Example: Spherical harmonics

Anisotropic line width

prm **p1** 0.01 min 0.0001

spherical_harmonics_hkl **sh1**

sh_order 2

lor_fwhm = **sh1 p1** / Cos(Th);

■ Anisotropic line asymmetry

prm **p2** 0.01 min 0.0001

spherical_harmonics_hkl **sh2**

sh_order 2


circles_conv = **sh2 p2** / Cos(Th);

p1, sh1, p2, sh2: Refineable parameters, arbitrary parameter names!

Anisotropic Refinement models



Example:
Set C:\Topas 4.2\tutorial\anisoLS\clay.inp in Launch mode
Run the Refinement
Edit the text file in text editor, uncomment Spherical_harmonics_hkl
lines and run again



spherical_harmonics_hkl sh
sh_order 8
 $\text{exp_conv_const} = (\text{sh}-1) / \sin(\text{Th});$

The Rietveld method

- is a standard-less full-profile approach to quantitative phase analysis
- uses every datapoint as a unique observation and least square methods to minimize the difference between calculated and measured intensities
- Residual of Least Square Refinement
$$R = \sum w_i(y_i - y_{c_i})^2$$

Rietveld Analysis requires:

- The crystal structure data for every phase in a mixture (unit cell and atomic positions)
- a model for the peak shapes and widths and a model for any aberrations
- a model for the background

The relative masses of all phases contributing to the diffraction pattern can be derived from the refinement using the simple relationship:

$$W_r = S_r (ZMV)_r / \sum_t S_t (ZMV)_t$$

W_r is the *relative* weight fraction of phase r in a mixture of t phases

S is the scale factor derived from Rietveld refinement

Z is the number of formula units per unit cell

M is the mass of the formula unit (atomic mass units)

V is the volume of the unit cell (\AA^3).

Measurement parameters:

Depend on sample properties: Determine peak width with a short test scan

- Adjust step size for sufficient amount of data points (6 data points above FWHM)
- Adjust time/step for sufficient counting statistics (at least a few thousand counts on the most intense peaks for basic quantitative analysis)

Typical measurement parameters for analysis of mineral samples:


- With Point detector:
 - 1° divergence slit
 - 4° Soller slits
 - 0.02 step size and at least 1 sec/step
- With 1-D detector:
 - 0.3° divergence slit and anti-airscatter screen
 - 4° Soller slits
 - 0.015 step size and at least 0.1 sec/step




- Good mixing, powder sample should be representative for bulk sample (no segregation)
- Smooth flat surface
- Sufficient particles to cover the x-ray beam area
- Sample thickness should be sufficient for correct intensities (infinite thickness for correct intensities)
- For low density materials thin layers will result in more accurate peak positions (transparency effect), often better for indexing
- Sufficient number of particles for good statistics (increase number of measured particles by rotation)

Sample Preparation


Powder Diffraction Basics



Diffraction of an ideal powder



Diffraction of textured materials



Diffraction of a small number of crystallites ("spotiness effect")

Sample Preparation

Powder Diffraction Basics

- Accurate, reproducible diffraction intensities require small crystallite size

typical intensity reproducibility for Quartz (113) reflection with $\text{CuK}\alpha$: is

15-20 μm	5-50 μm	5-15 μm	<5 μm
18.2%	10.1%	2.1%	1.2%

- The number of crystallites diffracting is related to size

diameter	40 μm	10 μm	1 μm
crystallites / 20mm^3	597.000	38.000.000	3.820.000.000
number diffracting	12	760	38.000

Smith, 1992

General Sample preparation requirements

Sufficient Grinding


■ Powder preparation

- Reduce particle size with mortar and pestle (or automatic mill) to smaller than 10 μm to increase number of crystallites and avoid spotiness effect
- Wet grinding with alcohol in MrCrone Mill is ideal for many geolocial samples.
- Too coarse grinding also effects preferred orientation




General requirements

Sufficient Grinding



Comparison of X-ray powder-diffraction patterns of a sample **mechanically ground under ethanol to <5 μm (blue)** and the same sample **dry ground by hand with a mortar and pestle (red)**. Note the large differences in intensity.

Rietveld Refinement Corrections



- Highlight corrections
- Use either Zero Error or Sample displacement as a peak shift correction
- Always Check Lorentz- Polarization correction
 - Type in the Bragg angle of the monochromator
 - Type 0, if no monochromator is present (Ni-filter, Sol-x detector)
- Depending on the sample you may refine Absorption correction

Using the correct structure file

Phase identification



- Correct phase identification is often more difficult than performing the actual Rietveld refinement.
- Use EVA and any complementary data that are available (e.g. XRF, EDX, thermal analysis, etc) to correctly identify all phases.
- Try to find the corresponding structure file in a database, not all structures are known.
- For complex samples with clays, a clay fraction may be necessary to correctly identify and quantify all phases in the sample.
- Quantitative Rietveld results are always normalized to 100% → If any phases are missed or mis-identified, it will lead to errors in the other phases as well.
- If amorphous phases are present, they can sometimes only be detected by using an internal standard.

Using the correct structure file

Structure databases






- Inorganic crystal structure database (ICSD) available at www.nist.gov or at <http://icsdweb.fiz-karlsruhe.de/index.php>
- Bruker-AXS structure database (a few hundred minerals)

Free databases:


- www.crystallography.net
- American Mineralogist Crystal structure database
<http://rruff.geo.arizona.edu/AMS/amcsd.php>
- Min-Chryst
<http://database.iem.ac.ru/mincryst/>
- Search for the structure in Google

Some structures in free data bases may have to be reformatted, if they do not conform to the standard CIF format.



Structure files

Entering the structure into TOPAS



- Highlight the scan file in the parameter window (the window below or in the context menu (right mouse click) allow several options for entering structures.
- Add Structure: Enter a structure manually (lattice parameters and atomic positions). This is sometimes required if the structure is only known from a journal.
- Load STR(s): Load a structure file from the Bruker database
- Load CIF(s): The Crystallographic Interchange Format is a standard format for structures. Most databases or publications have this format for crystal structures. Note that not every structure file in cif format that is found online, loads without problems. Occasionally, cif files have to be modified in a text editor to load correctly
- Load INP: Load a structure in Input file format (typically from TOPAS refinements)

Structure files

Example: Chlorite

Unit cell parameters



TOPAS - C:\users\old data\6008 RC\RC3-3_co radiation_3.pro - [RC 3-3 co rad 35KV 40ma 0.7div 0.2det sol-x.raw_1]

File View Fit Launch Tools Window Help

Structure Peak Type hkl Additional Convolutions Rpt/Text

	Use	Value	Code	Error	Min	Max
Use Phase	<input checked="" type="checkbox"/>					
Spacegroup	<input checked="" type="checkbox"/>	12				
a (Å)	<input checked="" type="checkbox"/>	5.3442000	a_chlorite	0.0000000	5.1900	5.5100
b (Å)	<input checked="" type="checkbox"/>	9.2660321	b_chlorite	0.0000000	8.9890	9.5450
c (Å)	<input checked="" type="checkbox"/>	14.2592097	c_chlorite	0.0000000	13.8420	14.6980
beta (°)	<input checked="" type="checkbox"/>	97.04304	be_chlorite	0		
Scale	<input checked="" type="checkbox"/>	2.08979475	sc_chlorite	0		
Cry Size						
Cry size L (nm)	<input checked="" type="checkbox"/>	360.2	cs_chlorite	0.0	32	
Cry size G (nm)	<input type="checkbox"/>	200.0	Refine	0.0		
LVol-IB (nm)	<input type="checkbox"/>	0.000		0.000	k:	1
LVol-FWHM (nm)	<input type="checkbox"/>	0.000		0.000	k:	0.89
Strain						
Strain L	<input type="checkbox"/>	0.1	Refine	0		
Strain G	<input type="checkbox"/>	0.1	Refine	0		
e0	<input type="checkbox"/>	0.00000		0.00000		
Wt% Rietveld	<input type="checkbox"/>	4.883		0.000		
Wt% of Spiked	<input type="checkbox"/>	0.000				
Wt% in Spiked sample	<input type="checkbox"/>	2.361		0.000		
Wt% in Original sample	<input type="checkbox"/>	2.574		0.000		
Cell Mass	<input type="checkbox"/>	1191.935		0.000		
Cell Volume (Å ³)	<input type="checkbox"/>	700.78128		0.00000		
Cry LAC (1/cm)	<input type="checkbox"/>	140.716		0.000		
Cry Density (g/cm ³)	<input type="checkbox"/>	2.824		0.000		
R Bragg	<input type="checkbox"/>	4.985				

Save Structure in STR format
View/Hide Structure
Create hkl's phase
Delete Structure
Paste INP to Node/Selections

- Values in red are refinable parameters (double click on the code column will change the state to „Fixed“)
- Constraining the lattice parameters and crystallite sizes to reasonable values helps with the refinement, especially for impurity phases or phases with broad peaks.
- Average crystallite size values for most minerals are typically above 30 nm and below 1000nm. Clays phases may require smaller crystallite sizes to fit broad peaks.
- Entering the weight of an internal standard in % allows quantification of amorphous materials

Structure files

Example: Chlorite

Atomic sites



TOPAS - C:\users\old data\6008 RC\RC3-3_co radiation_3.pro - [RC 3-3 co rad 35KV 40ma 0.7div 0.2det sol-]

File View Fit Launch Tools Window Help

Global


- RC 3-3 co rad 35KV 40ma 0.7c
 - Emission Profile
 - Background
 - Instrument
 - Corrections
 - Miscellaneous
 - Structures/ hkl Phases
 - Zincite
 - Quartz
 - Hematite
 - Magnetite
 - Muscovite 2M1
 - Chlorite IIb
 - Sites
 - Preferred Orientation
 - Str Output

Add Site(s) before selected site(s)
Add Site at bottom
Add Atom at selected site(s)
Paste INP to Node/Selections

Site	Np	x	y	z	Atom	Occ.	Beq.
1	s1	2	0.00000	0.00000	MG+2	0.3358	1
2					FE+2	0.6642	1
3	s2	4	0.00000	0.66780	MG+2	0.9703	1
4					FE+2	0.02971	1
5	s3	4	0.31500	0.00000	O-2	1	1
6	s4	8	0.18900	0.16670	O-2	1	1
7	s5	8	0.22480	0.16690	SI+4	0.66	1
8					AL+3	0.34	1
9	s6	4	0.80300	0.00000	O-2	1	1
10	s7	8	0.51100	0.22800	O-2	1	1
11	s8	4	0.82800	0.00000	O-2	1	1
12	s9	8	0.13100	0.34630	O-2	1	1
13	s10	4	0.00000	0.83360	MG+2	0.7016	1

Site	Np	x	y	z	Atom	Occ.	Beq.
1	s1	2	=0	=0	MG+2	=(1-p2)	Fix
2					FE+2	p2	Fix
3	s2	4	=0	!y2_chlori=0	MG+2	=(1-fe)	Fix
4					FE+2	fe	Fix
5	s3	4	!x3_chlori=0	!z3_chlori	O-2	Fix	Fix
6	s4	8	!x4_chlori!y4_chlori!z4_chlori	O-2	Fix	Fix	Fix
7	s5	8	!x5_chlori!y5_chlori!z5_chlori	SI+4	Fix	Fix	Fix
8					AL+3	Fix	Fix
9	s6	4	!x6_chlori=0	!z6_chlori	O-2	Fix	Fix
10	s7	8	!x7_chlori!y7_chlori!z7_chlori	O-2	Fix	Fix	Fix
11	s8	4	!x8_chlori=0	!z8_chlori	O-2	Fix	Fix
12	s9	8	!x9_chlori!y9_chlori!z9_chlori	O-2	Fix	Fix	Fix
13	s10	4	=0	!y10_chlo=1/2	MG+2	=(1-s)	Fix

- For quantitative mineral analysis the atomic sites (x,y,z) of the phase are never refined.
- In solution phases such as Chlorite the occupation factors for the Mg/Fe substitution may be refined.
- The x, y or z positions in blue color designate special positions for this unit cell.
- The atom or ion must be part of the drop down list in the atom column (with imported cif files this could lead to an error message)
- Parameter names like „p“ or „s“ have to be unique within a TOPAS project. If not, it will lead to this error message → rename parameters





Quantitative Rietveld Analysis

Exercise: Round-Robin example

File: C:\Topas 4.1\tutorial\qparr\cpd2.raw

Measurement circle radius: 173mm
Emission profile: CuKα5.lam
Background order 3, 1/x function
Divergence: 1°, Receiving slit: 0.3mm
Soller slits: prim. 4.6°, sec. 4.6°
LP factor: 26.6 (graphite monochromator)

Corundum, Brucite, Fluorite and Zincite Structures in str format in same directory


	hkl	Use	Value	Code
Direction 1	0 0 1	<input checked="" type="checkbox"/>	1	Refine
Direction 2		<input type="checkbox"/>	1	Refine
Fraction Dir. 1			0.5	Refine

Highlight Preferred Orientation in Brucite structure

Toggle Code column to Refine and Run Refinement again

Quantitative Rietveld Analysis


Exercise: Round-Robin example



Refined Lattice parameters and crystallite size for each phase

You may try to further improve the fit by introducing strain

Residuals and Goodness of Fit



Save your results in a TOPAS project file (.pro)

Quantitative Rietveld Analysis

More Examples:

Cement Clinker

File: C:\Topas 4.1\tutorial\qa\8486.raw

Measurement circle radius: 217.5mm

Emission profile: CuKα2.lam

Divergence: 1° , Receiving slit: 0.2mm

Soller slits: prim. 4° , sec. 4°

LP factor: 0 (Nickel filter)

Structures: C3Smmono_nishi.str, C3acub.str,
C2sbmu.str, C4af.str, CaO.str, Periclase.str

Iron Oxides

File: C:\Topas 4.1\tutorial\examples\Iron oxides
iron oxides.raw

Measurement circle radius: 217.5mm

Emission profile: CuKα2.lam

Divergence: 1° , Receiving slit: 0.2mm

Soller slits: prim. 4° , sec. 4°

LP factor: 0 (Sol-x detector)

Structures: in same directory

Aggregate with 11.4% Zincite as internal standard

File: C:\Topas 4.1\tutorial\Examples\Aggregate\aggregate with 11.4ZnO.raw

Measurement circle radius: 217.5mm

Emission profile: CuKα2.lam

Divergence: 0.3° , LynxEye detector with 3° opening

Soller slits: prim. 4° , sec. 4°

Structures: determine phases with Eva and use structures in c:\topas 4.1\structure database\

Natural Granodiorite

File: C:\Topas 4.1\tutorial\examples \cpd-granodiorite\granodio.raw

Measurement circle radius: 173mm

Emission profile: CuKα2.lam

Divergence: 1° , Receiving slit: 0.3mm

Soller slits: prim. 4.6° , sec. 4.6°

LP factor: 26.6 (graphite)

Structures: in same directory, Quartz, Feldspar, Albite, Biotite, less Chlinochlore, Hornblende, Zircon

Quantitative Rietveld Analysis

More Examples:

Paint panel

File: C:\Topas 4.1\tutorial\Examples\Paint\paint panel.raw

Measurement circle radius: 217.5mm

Emission profile: CuK α 2.lam

Divergence: 0.3° , LynxEye detector with 3° opening

Soller slits: prim. 4° , sec. 4°

Structures: determine phases with Eva and use structures in c:\topas 4.1\structure database\

Shale sample

File: C:\Topas 4.1\tutorial\Examples\Shale\Shale D4 06mm div 4dg soller 3dg LE_10pc ZnO.raw

Measurement circle radius: 201.5mm

Emission profile: CuK α 2.lam

Divergence: 0.3° , LynxEye detector with 3° opening

Soller slits: prim. 4° , sec. 4°

Structures: determine phases with Eva and use structures in c:\topas 4.1\structure database\


General error sources and trouble shooting:

- Check crystallite size of all phases for reasonable values. (it should not converge to zero or infinity), Constrain with min and max values.
- Check visually by highlighting the individual phase names, if trace phases are really present.
- Avoid correcting for preferred orientation when peaks of different phases are overlapping and both phases have preferred orientation (always try to minimize preferred orientation by better sample preparation)
- For complex multiphase sample, release parameters subsequently.
- Use internal standard to refine for displacement error, then fix it and only refine lattice parameters of other phases.

Corrections for non-ideal powder samples

Texture or preferred orientation

Slate Tile, as an example for strong texture




- Preferred orientation of needle or plate-like crystallites can often be corrected with the [March-Dollase](#) function
- For minerals with multiple cleavage directions a [spherical harmonics correction](#) can be applied.
- Acceptable limits of preferred orientation corrections should ideally be tested with a standard mixture

Illite, Kaolinite, Quartz


Correction for non-ideal powder samples

Texture or preferred orientation



← Unground sample

- Comparison of slate tile and ground slate (Calculated Illite patterns are highlighted)
- The (001) directions of Illite, Kaolinite and Paragonite were corrected with the March-Dollase model for both sample




← Ground sample

- Preferred Orientation is too severe to be corrected with March Dollase model
- Samples measured with Vantec-1 detector
0.1sec/step

Preferred orientation correction



Spherical harmonics correction using too many parameters





Negative peak intensities as a result of using the spherical harmonics correction for preferred orientation for Kaolinite

Corrections for non-ideal powder samples

Large grains lead to spotiness effect and Micro absorption
(demonstrated with area detector and weighted test mixture)



- Coarse powder of Albite, Microcline and Quartz (some grains are larger than 50 μ m)




- Identical Powder, wet ground in Micronising Mill for 30 minutes

Gadds data at 25cm sample-detector distance

Corrections for non-ideal powder samples

Large grains, spottiness effect, Micro absorption




← Unground sample

Weighted compositions:

Albite	29.64 wt. %
Microcline	36.63 wt. %
Quartz	33.73 wt. %

- For both samples, preferred orientation for Albite and Microcline were refined



← Ground sample

Gadds data at 25cm
sample-detector distance


Corrections for non-ideal powder samples

Micro absorption

Brindley correction



Unground sample with Brindley correction



Weighted compositions:

Albite	29.64 wt. %
Microcline	36.63 wt. %
Quartz	33.73 wt. %

■ Larger Quartz grains of 50µm were assumed.

- Instrument
- Corrections
- Miscellaneous
- Structures/ hkl Phases
 - Quantitative
 - SG-Lattice Parameters
 - Brindley Correction (Spherical Particles)
 - Display
- Quartz
 - Sites
 - Preferred Orientation
 - Str Output
- Microcline maximum

Brindley Correction (Spherical Particles) | Rpt/Text

Phase Name	Use	Radius(cm)	Pk Den.	Mass Ab.Coeff(cm ² /g)
Quartz	<input checked="" type="checkbox"/>	0.05	0.4	36.01106
Microcline maximum	<input checked="" type="checkbox"/>	0.005	0.4	49.87407
Albite intermediate	<input checked="" type="checkbox"/>	0.005	0.4	33.9074

Brindley Correction (Spherical Particles)

Phase Name	Use	Radius(cm)	Pk Den.	Mass Ab.Coeff(cm ² /g)
Quartz	<input checked="" type="checkbox"/>	0.05	0.4	36.01106
Microcline maximum	<input checked="" type="checkbox"/>	0.005	0.4	49.87407
Albite intermediate	<input checked="" type="checkbox"/>	0.005	0.4	33.9074

Note that the error would be considerably larger for phases with higher absorption contrast!

Corrections for non-ideal powder samples




- Preferred orientation effects
 - March Dollase
 - Spherical harmonics
- Micro absorption
 - Brindley correction for spherical particles
- Surface roughness
 - Correction according to Pitschke or Suortti

Corrections for these effects are available in TOPAS, but it is best to avoid or minimize these effects by good sample preparation!

- Flat samples with powder grain sizes below approx. 5 μ m.
- If possible, minimize preferred orientation by using back or side loaded sample holders

Quantifying Amorphous Phases Traditional Rietveld Method



Quantifying Amorphous Phases

Traditional Rietveld Method

Benefits

- Requires no standards or calibration
- More than one amorphous phase can be analyzed

Limitations

- Direct method
- Cannot correct for microabsorption errors
- Some amorphous material will not have a representative crystal structure
 - Available crystal structures (with long-range order) may not accurately represent material which only has short-range order (e.g. glasses).

Quantifying Amorphous Phases

Internal Standard Method

- The sample is "spiked" with a known mass of standard material and the QPA normalized accordingly
- The weight fractions of the crystalline phases present in each sample are estimated using the Rietveld methodology
- Concentrations to be corrected proportionately according to:

$$Corr(W_{\alpha}) = W_{\alpha} \frac{STD_{known}}{STD_{measured}}$$

where $Corr(W_{\alpha})$ is the corrected weight percent, STD_{known} the weighed concentration of the standard in the sample and $STD_{measured}$ the analyzed concentration

- The amount of amorphous material $W_{amorphous}$ can then be derived from:


$$W_{amorphous} = 1 - \sum_{j=1}^n Corr(W_j)$$

Quantifying Amorphous Phases

Internal Standard Method



Determination of amorphous component with internal standard



Phase	Wt.% in spiked sample	Wt.% in original sample	Wt.% det. without standard
Diopside	4.98	5.52	4.82
Forsterite	6.40	7.09	5.43
Anorthite	7.18	7.95	9.11
Volcanic glass	71.76	79.45	80.64

Quantifying Amorphous Phases

Internal Standard Method

Benefits

- Indirect method
- The Internal Standard Method is enabled in many Rietveld analysis packages

Limitations

- Only the sum of all amorphous and unidentified phases can be reported
- Cannot correct for microabsorption errors
- The sample is contaminated
- The standard addition process is laborious (weighing, mixing), and not feasible in industrial, automated sample preparation environments
- The method relies upon obtaining a standard of appropriate absorption contrast to prevent the introduction of a microabsorption problem

Quantifying Amorphous Phases

PONKCS Method

- Phases with Partial Or No Known Crystal Structure are characterized by measured rather than calculated structure factors
- Follows the same general form as that used in the Rietveld Method but now includes all crystalline and amorphous phases characterized by either calculated or empirical structure factors
- For all phases α using empirically derived structure factors ZMV "calibration constants" must be derived, e.g. via an internal standard s

$$(ZMV)_{\alpha} = \frac{W_{\alpha}}{W_s} \frac{S_s}{S_{\alpha}} (ZMV)_s$$

- A one time calibration per phase with a single standard mixture is usually sufficient

Quantifying Amorphous Phases

PONKCS Method

Benefits


- The amorphous phase is included in the analysis model
- More than one amorphous phase can be analyzed
- Potential to minimize errors related to microabsorption, when ZMV's have been calibrated for all phases
 - Consider to use the PONKCS methodology also for all crystalline phases with known crystal structures are known
 - It is highly recommended to create a database with "PONKCS phases"; such a database can be used in full analogy to crystal structure databases

Limitations

- Direct method
- Requires availability of a standard mixture to derive an empirical ZMV

Quantifying Amorphous Phases


PONKCS Method-walk through for Compritol



- Scan of pure Compritol
- Step 1: try to find a "unit cell" that will fit the scan. If not enough peaks are present for LSI algorithm, try LP search
- Step 2: after getting a satisfactory fit, fix all peak intensities and uncheck "Delete hkl on refinement"

Quantifying Amorphous Phases


PONKCS Method-walk through for Compritol



- Scan of pure Compritol
- Step 1: try to find a "unit cell" that will fit the scan. If not enough peaks are present for LSI algorithm, try LP search
- Step 2: after getting a satisfactory fit, fix all peak intensities and uncheck "Delete hkl on refinement" and LP search

Quantifying Amorphous Phases


PONKCS Method-walk through for Compritol



- Scan of 50% mix of Compritol and known standard (Corundum)
- Step 3: replace scan of pure compound with known mixture
- Step 4: Constrain all lattice parameters and crystallite size values with reasonable limits and refine scaling factor

Quantifying Amorphous Phases


PONKCS Method-walk through for Compritol



- Scan of 50% mix of Compritol and known standard (Corundum)
- Step 3: Empirically determine a value for "mass" to match the known concentration of the mixture

Quantifying Amorphous Phases


PONKCS Method-walk through for Compritol



- Step 4: Once that value is known, save the hkl_Compritol phase (in .inp format). It can now be used like a normal structure in any other mixture

Quantifying Amorphous Phases


PONKCS Method-walk through for Prosolv



- Scan of pure Prosolv
- Step 1: try to find a "unit cell" that will fit the scan. If not enough peaks are present for LSI algorithm, try LP search
- Step 2: after getting a satisfactory fit, fix all peak intensities and uncheck "Delete hkl on refinement"

Quantifying Amorphous Phases

PONKCS Method-walk through for Prosolv




- Scan of pure Prosolv
- Step 1: try to find a "unit cell" that will fit the scan. If not enough peaks are present for LSI algorithm, try LP search
- Step 2: after getting a satisfactory fit, fix all peak intensities and uncheck "Delete hkl on refinement" and LP search

Quantifying Amorphous Phases


PONKCS Method-walk through for Prosolv

- Scan of 50% mix of Prosolv and known standard (Corundum)
- Step 3: replace scan of pure compound with known mixture



Quantifying Amorphous Phases


PONKCS Method-walk through for Prosolv



- Scan of 49.1% mix of Prosolv and known standard (Corundum)
- Step 3: Once that empirical value is known, save the hkl_Compritol phase (in .inp format). It can now be used like a normal structure in any other mixture

Quantifying Amorphous Phases


PONKCS Method-walk through for Klucel



- Scan of pure Klucel
- Step 1: try to find a "unit cell" that will fit the scan. If not enough peaks are present for LSI algorithm, try LP search
- Step 2: after getting a satisfactory fit, fix all peak intensities and uncheck "Delete hkl's on refinement"

Quantifying Amorphous Phases


PONKCS Method-walk through for Klucel



- Scan of pure Klucel
- Step 1: try to find a "unit cell" that will fit the scan. If not enough peaks are present for LSI algorithm, try LP search
- Step 2: after getting a satisfactory fit, fix all peak intensities and uncheck "Delete hkl on refinement" and LP search

Quantifying Amorphous Phases


PONKCS Method-walk through for Klucel



- Scan of 50% mix of Klucel and known standard (Corundum)
- Step 3: replace scan of pure compound with known mixture


Quantifying Amorphous Phases

PONKCS Method-walk through for Klucel



- Scan of 60% mix of Klucel and known standard (Corundum)
- Step 3: Once that empirical value is known, save the hkl_klucel phase (in .inp format). It can now be used like a normal structure in any other mixture

Quantifying Amorphous Phases PONKCS Method-Test mixture




- Test mixture
 - 39% Ibuprofen
 - 18.7% Prosolv
 - 18% Compritol
 - 24.24% Klucel
- Background was determined on pure Ibuprofen sample to avoid correlation effects with broad excipient peaks
- Ibuprofen structure was also refined on pure sample to minimize preferred orientation effects → no further correction necessary

Industrial application: Cement + flyash

Rietveld refinement with TOPAS, only crystalline components



Sample 7




- LynxEye Detector
- 0.3° divergence, 4° Soller slits
- 0.02° /step, 0.2sec/step
- Measurement time: 9 minutes

A high order polynomial function was necessary to fit the background, which is an indication for amorphous material in the sample.

Industrial application: Cement and flyash

Determination of amorphous component

Method 1: spiking the sample with an internal standard



Background was fitted with high-order polynomial function

8.27 wt% Zincite was added to the sample as an internal standard

The elevated background is modeled with high-order polynomial function.

The sample is spiked with a known amount of an internal standard.

Structure	Peak Type	hkl	Additional Convolution
sample 7_500mm 0.3div			
Emission Profile			
Background			
Instrument			
Corrections			
Miscellaneous			
Structures/ hkl Phases			
C3S monoclinic (NISHI)			
C2S beta (MUMME)			
C3A cubic			
C4AF			
Quartz			
Gypsum			
Periclase			
Zincite			
Sites			
Preferred Orientation			

Structure	Peak Type	hkl	Additional Convolution
Use Phase	<input checked="" type="checkbox"/>		
Spacegroup		186	
a (Å)		3.2495907	a_zincite
c (Å)		5.2067137	c_zincite
Scale	<input checked="" type="checkbox"/>	0.00369548	sc_zincite
Cry Size			
Cry Size L(nm)	<input checked="" type="checkbox"/>	153.5	cs_zincite
Cry Size G(nm)	<input type="checkbox"/>	200.0	Refine
LVolHB(nm)	<input type="checkbox"/>	0.000	
LVolFWHM(nm)	<input type="checkbox"/>	0.000	
Strain			
Strain L	<input type="checkbox"/>	0.1	Refine
Strain G	<input type="checkbox"/>	0.1	Refine
σ0	<input type="checkbox"/>	0.00000	
Wt% Rietveld		13.951	
Wt% of Spiked	<input checked="" type="checkbox"/>	8.270	

Industrial application: Cement + flyash

Determination of amorphous component

Method 1: spiking the sample with an internal standard



sample 7_500mm 0.3div

- Emission Profile
- Background
- Instrument
- Corrections
- Miscellaneous
- Structures/ hkl Phases
 - Quantitative**
 - SG-Lattice Param
 - Brindley Correction
 - Display
- C3S monoclinic (NIST)
- C2S beta (MUMME)
- C3A cubic
- C4AF
- Quartz
- Gypsum
- Periclase

Paste INP to Node/Selections

Phase Name	Wt% Rietve	Use	Wt% of Spil	Wt% in Spike	Wt% in Origina	Cell Mass	Cell Vol (Å ³)
C3S monoclinic	55.459	<input type="checkbox"/>	0.000	32.875	35.839	8219.355	4307.46073
C2S beta (MUMME)	12.849	<input type="checkbox"/>	0.000	7.617	8.304	688.953	347.61380
C3A cubic	5.123	<input type="checkbox"/>	0.000	3.037	3.310	6484.622	3542.47539
C4AF	9.019	<input type="checkbox"/>	0.000	5.346	5.828	971.922	429.80950
Quartz	0.000	<input type="checkbox"/>	0.000	0.000	0.000	180.252	112.81456
Gypsum	3.599	<input type="checkbox"/>	0.000	2.134	2.326	672.560	495.62576
Periclase	0.292	<input type="checkbox"/>	0.000	0.000	0.000	161.217	74.79168
Zincite	13.951	<input checked="" type="checkbox"/>	8.270	8.270	0.000	162.778	47.61586

Amorphous Content


	Value	Error
Wt% in Spiked sample	40.722	0.772
Wt% in Original sample	44.393	0.841

Phase amounts and amorphous content in original sample is calculated in the TOPAS software

Industrial application: Cement + flyash

Determination of amorphous component

Method 2: Modeling the background with „amorphous structure“




Background from
0-order
polynomial
function

Background from
amorphous
component

The elevated background from the amorphous component is not modeled with a polynomial but mostly with an amorphous “structure”. This was achieved by using a normal quartz structure with very small crystallite size and adjusting the lattice parameters to fit the amorphous halo.

This standard less method works reliable and can be automated easier but should be checked for each sample type with the internal standard method.



www.bruker-axs.com